



SKR058-MT

ANESTHESIA TROLLEY



SKSurgitech™
by Saikang Medical

Care you precisely

SKR058-MT

FEATURES AND HIGHLIGHTS

Material

The ABS material is an innovative engineering plastic. Antibacterial, moisture-proof, anti-corrosion, wear-resistant, non-toxic and tasteless, never rust. Washable, easy cleaning.



Design

The overall design and production are exquisite, the styling and color matching are reasonable, installation is simple, and the structure is stable and durable.



Drawer

Five layers of drawers, a variety of storage, more convenient and quick access.



Drawer compartment

The interior of the drawer is equipped with transparent acrylic classified medicine grids, which can be used for classifying and placing first aid medicines. The medicine grids can be freely combined.

Technical parameters

▪ Overall dimensions	820W x 640D x 1080H mm
▪ Anesthesia rack	size600*300*320mm
▪ Caster	5"two casters with brake, two without brake
▪ Work board	340*250mm
▪ Packing size	820x650x1150 mm
▪ Productweight	65 kg

Technical configuration

- 8 lateral medicine hanging boxes
- Drawer divider system
- 20L trash can with lid
- Sealed handle

Accessories

Anesthesia rack

Double-row adjustable eight-hook anesthesia stand with a sturdy and durable aluminum alloy support and 8 acrylic material drug hanging boxes.

Drawer slide

The drawer uses high-strength three-section silent damping slides. When the drawer is closed to the final distance, hydraulic pressure is used to slow down the speed and reduce the impact force.

Working table

The trolley is equipped with a collapsible working table, which can be used for writing; Work board size: 340*250mm;

Sealed handle

The trolley is equipped with a sealed handle, which can lock five drawers at the same time.

Trash can

Clamshell ABS trash can on the right side, specification: 20L

Caster

At the bottom, it uses two 5-inch silent brake, anti-winding, high-end, silent casters, allowing for flexible movement and a high load-bearing capacity.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue