



More production and solidarity against Covid-19

Quality test laboratory	4
Sundew Day Surgery	5
Innovating: D8d bed	6

We increase our production against the pandemic

Electric and manual beds, stretchers, complementary furniture, ventilators, masks... all our products are being of great help throughout the world in this difficult time.

Any professional activity and selfless help is very necessary at this time that we are living in a pandemic. In these months we have increased our production and sales, especially in products and materials related to the hospital fight against Covid-19: electric and manual beds, stretchers, complementary furniture, ventilators, masks, etc. Thanks to our large production capacity and our stock in stock, our customers and users are receiving the material as urgently as possible and putting it into use quickly, since time counts now more than ever.

As a result, if we already regularly receive evaluations and images of our products from many countries, now all this has also multiplied. And there is nothing more rewarding than seeing all expectations met. If we are already very satisfied with the work done and with being able

to contribute to improving medical services around the world, at this delicate moment our pride is even greater.

Donated face masks

On the other hand, at Saikang Medical we also try to provide our help against the pandemic outside of our productive activity. We have donated almost 30,000 masks, distributing them among our clients in more than 20 countries. Thus, these weeks we have also received messages of thanks and some photos of our clients from various countries (Italy, Poland, Russia, Kenya, Bangladesh, Singapore, Thailand, Cambodia, Myanmar, the Philippines, etc.).

Our most sincere thanks to all of you. We will continue to contribute in everything that is in our hands. Thank you all very much for helping us to export health!



The images sent by our clients in different countries show our products at their destination, most of them already installed and ready to use. In the last photos, some clients show us and appreciate the face masks donated by our company.



Different symptoms common cold, flu and Novel Coronavirus.



Symptoms caused by Common Cold:
Light symptoms without pneumonia

Symptoms caused by Flu:
Fever, cough and muscular pain

Symptoms caused by Novel Coronavirus:
Cough, weakness, chest distress, and difficult breathing



Why should anyone in close contact with COVID-19 infectors be under medical observation for 14 days?



The average incubation period is about 7 days. And the longest incubation is up to 12 days.
Then, Home-based medical observation is necessary.



When should we receive medical treatment?



Some symptoms like cough, chest distress, and difficult breathing appear.
And have a close contact with people from epidemic area.



What if you are diagnosed with a novel coronavirus?



Receive an isolation treatment

Provide private contact history



Wash hands with 7 procedures



We share some recommendations that can help us against the virus. We hope they are of interest to you and will help everyone.

New improvements in our quality testing laboratory

At Saikang Medical we continue to invest in equipment to complete our quality testing laboratory. Thus, we not only continue to guarantee that our certified products meet all standards, but also provide us with one of our production objectives, which is innovation for the continuous improvement of our products and components.

With the latest implementations, the experts in our quality department have a widening range of methods for control and anticipation: Webster and Rockwell hardness tests (to measure the hardness of aluminum, iron and steel and thus guarantee the quality of the metal parts of our products, such as the chassis, platform or guardrails of our hospital beds), air spring test, IV pole load capacity test, color difference test (whether for plastic, paint, ink printing or fabric), head and foot board tail stress test, 10° slope test (for the correct fastening of wheels and bed brakes on a sloping surface), coating thickness test, cross-cut test (to measure the resistance of paints on the bed surface), motor mounting seat test (to ensure the stability of a correct motor mount), electric bed function test, castor test (standard abrasion resistance for wheels under heavy load), threshold test (with obstacles to check mobility under loads), salt spray test (to measure the durability of metal parts) and vibration platform for packaging (to simulate transport conditions and ensure safe shipping of our products).



In the images, some of the tests that are carried out, such as the vibration platform for packaging, loads and weight mobility, strength and durability of all the components (such as the ends and wheels), etc.



Sundew Day Surgery, a new hospital in Australia

Saikang Medical has supplied complete equipment in Perts



This time we bring information from the continent of Oceania, where we have been able to contribute to equipping the new Sundew Day Surgery hospital, in Perth (Australia). The project has included 10 medical beds, 3 operating tables, 12 operating lights, and complementary medical equipment such as monitors, medical carts, office tables, and waiting chairs.



This health center belongs to the Nexus Hospitals group and is northern Perth's newest private day surgery. According to the heads of the center, "has been purpose built for short stay surgery, both day and overnight paediatric and adult patients, and provide a warm and welcoming environment with the latest technology and a highly specialised surgery staff".



As we see in the images received, the facilities are complete and are currently in operation. So, we are very happy to participate in this new project and satisfied with the work done, since the final evaluation of the client has been very positive.

D8d

The best bed option for ICU use in this pandemic period

The D8d model, belonging to the Grace Care series, is one of the models that is receiving the most requests in this complicated period. It is an electric multifunction bed that, in addition to covering all the needs for its perfect use for ICU (back-rest and knee-rest adjustment, reverse Trendelenburg and Trendelenburg, CPR, etc.), has a modern design and features that meet all the expectations of professionals.



PRODUCT VIDEO

Technical parameters

- Size (LxW): 2280x1090 mm
- Bed height adjustment: 435-725 mm
- Mattress platform: 1925x900 mm
- Safe working load: 250 kg
- Back-rest adjustment: 0-75° ($\pm 10^\circ$)
- Knee-rest adjustment: 0-35° ($\pm 10^\circ$)
- Reverse Trendelenburg: 0-12°
- Trendelenburg: 0-12°



Technical configuration

- Manual CPR
- One button stop
- One button trendelenburg
- Nurse lock function
- Linear motor: 4 pcs.
- Angle indicator: 2 pcs.
- 5" double side castors with central locking system: 4 pcs.
- Central locking pedal: 2 pcs.
- ABS side rail with integrated control panel: 1 set
- ABS bed ends with integrated LCD control panel: 1 set
- Wall bumper wheel ($\phi 85$ mm): 4 pcs.
- Trapeze socket: 2 pcs.
- IV pole: 1 set
- IV pole prevision: 4 pcs.
- Drainage hook: 2 pcs.

Optional

- Electrical CPR
- Weighing system with LCD screen
- Bed extension (20 cm)
- Back-up battery
- Cardiac position
- ABS bed surface
- X-ray translucent mattress platform
- Directional wheels function
- Monkey pole
- Aluminium balkan frame
- Basic foam mattress (12 cm thickness)
- Anti-decubitus foam mattress (14 cm thickness)
- Oxygen bottle holder
- Single arm lifting pole
- Utility shelf



In brief

Vitamin D seems to help against Covid-19

A new study has found an association between low vitamin D levels and high amounts of Covid-19 cases (and associated death rates), in 20 European countries. The research has been led by Dr. Lee Smith of Anglia Ruskin University and Dr. Petre Cristian Ilie of QEHL (Queen Elizabeth Hospital Kings Lynn NHS Trust), both entities in the UK. The human body can obtain vitamin D by exposure to sunlight, diet, or supplements. Previous observational studies noted an association between low vitamin D levels and susceptibility to acute respiratory tract infections. The SARS-CoV-2 coronavirus is known to cause excess proinflammatory cytokines.



New company video on Alibaba

Saikang Medical has been present for 14 years already on the Alibaba platform, the online giant of international trade business-to-business. From time to time we renew the promotional materials and last April a technician visited us to prepare a new presentation video for our company. With him we toured our facilities to record and be able to show the different areas and production processes in the best way.

Once it's ready, the new video will be visible in our two Alibaba stores, saikang.en.alibaba.com and saikang-surgical.en.alibaba.com. In addition, we will also share it on our social media.



Next appointment in Shanghai

As we know, most medical exposures have been and are being canceled or postponed throughout these months. Since Saikang Medical participated in Arab Health (Dubai) in late January, we have not been able to exhibit in any other city in the world.

Everything indicates that the next appointment will be at China International Medical Equipment, scheduled for October 19-22. This exhibition was already postponed in June, but we hope it can be held on this new date and we can meet at the National Exhibition and Convention Center in Shanghai. We will inform you as there are new confirmations.



SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

📞 +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue