

CELEBRATING
ANNIVERSARY

24th

· 2002-2026 ·



Healthier your life

Advanced Electric Dialysis Chair Solution

www.saikangmedical.com

SAIKANG Medical®

- Optimizing Dialysis Workflow
- Ensuring Safety
- Improving Care Quality

STRENGTH OF SAIKANG

24

24 years focusing on design, production, service and technical consultation of medical care equipment

150

Located in Zhangjiagang City, 150km nearby Shanghai Economic Center

1000000

Saikang delivered 1,000,000 sets of products globally. Products all over Asia, Middle East, North America and Africa

400

There are more than 400 employees worldwide, more than 50 qualified engineers, and more than 100 trained technicians

245

Equiped with 245 sets of various automatic production facilities, Saikang follow strict standard control system on every production process, enable the traceable of each part

50000

The finished phase I production area covers more than 30,000 square meters, and phase II will be more than 20,000 square meters of advanced production area

31

Saikang offer comprehensive medical solutions with 31 product categories

562

Safe, efficient and easy-to-use medical solution including emergency, nursing and ICU wards of more than 500 series of products

32

With CE certification, ISO9001/13485 quality system certification, U.S. FDA and other 32 certifications

189

Ownership of 189 patents, continuously leading new standards in defining product quality and service efficiency

150

Since the establishment of the company, products and solutions have been used in more than 150 countries

1000

Saikang products are committed to creating a modern nursing environment, accumulatively serving more than 1,000 medical and health care organizations



LEADING BRAND OF MEDICAL FURNITURE

Patients in medical rehabilitation often have to face huge challenges. Many people who suffering from diseases or accidents are faced with retraining their minds and bodies and must adapt to the new environment. In this case, the positive atmosphere would offer decisive help. Hospital beds and furniture should provide the best support for patients and nurses. Considering this case, Saikang devoted to provide medical care solutions and perfect choices for medical institutions and rehabilitation centers, to help patients achieve rapid or long-term recovery.

Saikang was founded in 2002, now is a leading global provider of medical furniture and solutions. Head office locates in Suzhou, and branches in US, Dubai, Russia, Brazil, Thailand and Kenya, having 400 employees worldwide. Saikang has formed a large global network of research and development, marketing and service. Currently, Saikang products and medical solutions have been used in more than 149 countries and regions in the world.

Common Challenges in Dialysis Care



Limited Support for Different Body Types

Non-ergonomic seating designs fail to support different body sizes and postures, leading to pressure points and poor sitting tolerance.



Safety Risks During Adjustment

Unstable structures or sudden movements during adjustment may increase the risk of patient discomfort or injury.



Long Treatment Duration

Dialysis treatments often last several hours, causing discomfort and fatigue when patients remain in the same position for extended periods.

Difficult Position Adjustment

Manual position adjustment requires physical effort from caregivers and may interrupt treatment flow, reducing efficiency and safety.

High Cleaning & Infection-Control Demands

Dialysis environments require frequent cleaning and strict infection control, making complex structures and hard-to-clean surfaces a challenge.



Why Choose Our Electric Dialysis Chair?

Multi-Position Electric Adjustment

Backrest, leg rest, and seat positions can be adjusted electrically to support frequent posture changes, improving comfort during long dialysis sessions.



Ergonomic Seating Design

Ergonomically contoured seating provides balanced support for different body types, reducing pressure and improving sitting tolerance.



Stable & Safe Structure

A reinforced frame and stable base ensure smooth and secure adjustment in all positions, enhancing patient safety throughout treatment.



Caregiver-Friendly Operation

Easy-to-use hand controls reduce physical effort for caregivers and help maintain an efficient treatment workflow.



Infection-Control Ready Design

Smooth, sealed surfaces allow fast cleaning and disinfection, meeting the strict infection-control standards of dialysis centers.



Integrated Accessories & Options

Optional accessories enable customization for different clinical requirements, improving overall treatment convenience.

Coordinated Dialysis Environment

Designed to integrate seamlessly with dialysis room furniture and equipment, creating a clean and organized treatment environment.



TECHNICAL SPECIFICATIONS



Model	SKE-131	SKE-132	SKE-133	SKE-134
Size (L×W)	1800×580±10mm	1750×580±10mm	2000×580±10mm	1930×580±10mm
Backrest	770×580±10mm	750×580±10mm	770×580±10mm	780×580±10mm
Seat	500×580±10mm	460×580±10mm	470×580±10mm	480×580±10mm
Leg	(460-580)×470±10mm	400×470±10mm	400×540±10mm	580×530±10mm
Footboard	380×220±10mm	/	430×210±10mm	240×380±10mm
Pillow	360×240±10mm	/	/	/
Seat Height	560±10mm	520±10mm	/	/
Mattress Thickness	80±10mm	80±10mm	80±10mm	80±10mm
Back-rest Adjustment	-10°~70°±5°	-8°~70°±5°	-8°~70°±5°	-8°~70°±5°
Knee-rest Adjustment	-50°~0°±5°	-70°~0°±5°	-60°~0°±5°	-50°~0°±5°
Armrests Swing Angle	0°~180°	0°~180°	0°~180°	0°~180°
Height Adjustment	/	/	550~700	550~700
Footboard Adjustable Stroke	/	/	/	150
Armrests Swing Angle	/	/	/	/
With CPR Function	Yes	No	Yes	Yes
Hand Controller	1pc	1pc	1pc	1pc
Power Line	1pc	1pc	1pc	1pc
3-inch Brake Wheel	4pcs	4pcs	4pcs	4pcs
Arm Support	2pcs	2pcs	2pcs	2pcs
Small Dining Table	/	1pc	1pc	1pc
USB Port	/	/	/	/
Side Guide Rails	/	/	/	/
Plug-in Headrest	/	/	/	/
Seamless Soft-padded Seat Cushion	/	/	/	/
Silent Castors	/	/	/	/
Standard Pillow	/	/	/	/



Model	SKE-135	SKE-136	SKE-188A
Size (L×W)	1750×550±10mm	2020×600±10mm	2040×700±10mm
Backrest	760×590±10mm	610×600±10mm	820×700±10mm
Seat	480×550±10mm	500×600±10mm	470×700±10mm
Leg	(460-580)×370±10mm	570×600±10mm	660×700±10mm
Footboard	380×220±10mm	260×640±10mm	680×500±10mm
Pillow	360×240±10mm	470×260±10mm	360×230±10mm
Seat Height	540±10mm	/	/
Mattress Thickness	80±10mm	90±10mm	80±10mm
Back-rest Adjustment	-10°~70°±5°	0°~70°±5°	0°~70°±5°
Knee-rest Adjustment	-70°~0°±5°	-50°~0°±5°	-40°~0°±5°
Armrests Swing Angle	0°~180°	0°~180°	/
Height Adjustment	/	620~870	580~670
Footboard Adjustable Stroke	/	150	150
Armrests Swing Angle	/	6°	/
With CPR Function	No	Yes	No
Hand Controller	1pc	1pc	1pc
Power Line	1pc	/	1pc
3-inch Brake Wheel	4pcs	4pcs	/
Arm Support	2pcs	2pcs	/
Small Dining Table	/	1pc	/
USB Port	/	1pc	/
Side Guide Rails	/	2pcs	1 set of Foldable Aluminum Guardrail
Plug-in Headrest	/	1pc (with a small pillow)	/
Seamless Soft-padded Seat Cushion	/	Available	/
Silent Castors	/	/	4pcs
Standard Pillow	/	/	1pc

SKE-131

FEATURES AND HIGHLIGHTS



Technical parameters

▪ Size(L×W):	1800×580mm(±10mm)
▪ Back:	770×580mm(±10mm)
▪ Seat:	500×580mm(±10mm)
▪ Leg:	(460-580)×470mm(±10mm)
▪ Foot board:	380×220mm(±10mm)
▪ Pillow:	360×240mm(±10mm)
▪ Seat height:	560mm (±10mm)
▪ Mattress thickness:	80mm(±10mm)

Technical configuration

☑ Hand Controller	1 pc
☑ Power Line	1 pc
☑ 3-inch brake wheel	4 pc
☑ Arm support	2 pc

Functions

▪ Back-rest Adjustment:	-10°- 70°(±5°)
▪ Knee-rest Adjustment:	-50°-0°(±5°)
▪ Armrests swing out:	0°-180°
▪ With CPR	



QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Steels
Baosteel from Fortune Global 500.



Plastic
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Welding
Panasonic Robotic ensure 360° full smooth wedding.



Electric Motor
Linear motor, with UL/ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.

Accessories

Chair Cover
Optional genuine leather or supebric cover, with artisanal making, comfortable hand feeling.

Pillow
Ergonomics design removable pillow.

Handrail
Upgraded armrest, can be right-left, front-back adjustable. Donor can control universal position. Armrest is R120° concave curve, avoid donor's arm sliding, ensure the donor's safety and comfortable.

Dustproof Base Cover
With dustproof base cover, protect motors and lines, easy to clean.

Ultrafiber leather color table



* Note: the above is the physical scan, there may be a certain color difference in kind.

SKE-132

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Size(LxW)	1750×580 mm (±10mm)
▪ Back	750×580 mm (±10mm)
▪ Seat	460×580 mm (±10mm)
▪ Leg	400×470 mm (±10mm)
▪ Seat height	520 mm (±10mm)
▪ Mattress thickness	80 mm (±10mm)

Functions

▪ Back-rest Adjustment	-8° - 70° (±5°)
▪ Knee-rest Adjustment	-70° - 0° (±5°)
▪ Armrests swing out	0° - 180°



Technical configuration

☑ Hand Controller	1pc
☑ Power Line	1pc
☑ 3-inch brake wheel	4pc
☑ Arm support	2pc
☑ Small Dining table	1pc



QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Steels
Baosteel from Fortune Global 500.



Plastic
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.




Welding
Panasonic Robotic ensure 360° full smooth wedding.



Electric Motor
Linear motor, with UL/ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.

Accessories



Chair Cover
Optional genuine leather or supebric cover, with artisanal making, comfortable hand feeling.




Table Board
Adjustable ABS tray table, with antiskid convex edge, tank for cup and cellphone. USB/5V optional. Very convenient for donor.



Handrail
Upgraded armrest, can be right-left, front-back adjustable. Donor can control universal position. Armrest is R120° concave curve, avoid donor's arm sliding, ensure the donor's safety and comfortable.



Dustproof Base Cover
With dustproof base cover, protect motors and lines, easy to clean.

Ultrafiber leather color table



* Note: the above is the physical scan, there may be a certain color difference in kind.

SKE-133


FEATURES AND HIGHLIGHTS

Technical parameters

▪ Size(LxW):	2000 x 580 mm(±10mm)
▪ Back:	770 x 580 mm(±10mm)
▪ Seat:	470 x 580 mm(±10mm)
▪ Leg:	400 x 540 mm(±10mm)
▪ Foot board:	430 x 210 mm(±10mm)
▪ Mattress thickness:	80 mm(±10mm)

Functions

▪ Back-rest Adjustment	-8° - 70° (±5°)
▪ Knee-rest Adjustment	-60° - 0° (±5°)
▪ Height adjustment:	550 - 700 mm(±10mm)
▪ With CPR	
▪ Armrests swing out:	0° - 180°



Technical configuration

☑ Hand Controller	1pc
☑ Power Line	1pc
☑ 3-inch brake wheel	4pc
☑ Arm support	2pc
☑ Small Dining table	1pc



QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Steels
Baosteel from Fortune Global 500.



Plastic
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.




Welding
Panasonic Robotic ensure 360° full smooth weddling.



Electric Motor
Linear motor, with UL/ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.

Accessories



Chair Cover
Optional genuine leather or supebric cover, with artisanal making, comfortable hand feeling.



Foot Board
Foot pedal connects with leg rest, adjustable, can automatic extension the chair lengthh 250mm, realize both foot step and leg support functions.



Handrail
Upgraded armrest, can be right-left, front-back adjustable. Donor can control universal position. Armrest is R120° concave curve, avoid donor's arm sliding, ensure the donor's safety and comfortable.




Table Board
Adjustable ABS tray table, with antiskid convex edge, tank for cup and cellphone. USB/5V optional. Very convenient for donor.

Ultrafiber leather color table



* Note: the above is the physical scan, there may be a certain color difference in kind.

SKE-134


FEATURES AND HIGHLIGHTS

Technical parameters

▪ Size(LxW)	1930 × 580 mm(±10mm)
▪ Back	780 × 580 mm(±10mm)
▪ Seat	480 × 580 mm(±10mm)
▪ Leg	580 × 530 mm(±10mm)
▪ Foot board	240 × 380 mm(±10mm)
▪ Mattress thickness	80 mm(±10mm)

Functions

▪ Back-rest Adjustment	-8° - 70°(±5°)
▪ Knee-rest Adjustment	-50° - 0°(±5°)
▪ Height adjustment	550 - 700mm(±10mm)
▪ Foot board adjustable distance	150 mm(±10mm)
▪ With CPR	
▪ Armrests swing out	0° - 180°



Technical configuration

☑ Hand Controller	1pc
☑ Power Line	1pc
☑ 3-inch brake wheel	4pc
☑ Arm support	2pc
☑ Small Dining table	1pc



QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Steels
Baosteel from Fortune Global 500.



Plastic
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.




Welding
Panasonic Robotic ensure 360° full smooth wedding.



Electric Motor
Linear motor, with UL/ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.

Accessories



Chair Cover
Optional genuine leather or supebric cover, with artisanal making, comfortable hand feeling.





Table Board
Adjustable ABS tray table, with antiskid convex edge, tank for cup and cellphone. USB/5V optional. Very convenient for donor.



Handrail
Upgraded armrest, can be right-left, front-back adjustable. Donor can control universal position. Armrest is R120° concave curve, avoid donor's arm sliding,ensure the donor's safety and comfortable.



Foot Board
150mm foot board adjustable distance.

Ultrafiber leather color table



* Note: the above is the physical scan, there may be a certain color difference in kind.

SKE-135

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Size(LxW):	1750mm×550mm(±10mm)
▪ Back:	760×590mm(±10mm)
▪ Seat:	480mm×550mm(±10mm)
▪ Leg:	(460-580)×370mm(±10mm)
▪ Foot board:	380×220mm(±10mm)
▪ Pillow:	360×240mm (±10mm)
▪ Seat height:	540mm (±10mm)
▪ Mattress thickness:	80mm(±10mm)
▪ Back-rest Adjustment:	-10°-70°(±5°)
▪ Knee-rest Adjustment:	-70°-0°(±5°)
▪ Armrests swing out:	0°-180°



Technical configuration

☑ Hand Controller:	1pc
☑ Power Line:	1pc
☑ 3-inch brake wheel:	4pc
☑ Arm support:	2pc



QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Steels
Baosteel from Fortune Global 500.



Plastic
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Welding
Panasonic Robotic ensure 360° full smooth welding.




Electric Motor
Linear motor, with UL/ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.

Accessories



Chair Cover
Optional genuine leather or supebric cover, with artisanal making, comfortable hand feeling.



Pillow
Ergonomics design removable pillow.



Handrail
Upgraded armrest, can be right-left, front-back adjustable. Donor can control universal position. Armrest is R120° concave curve, avoid donor's arm sliding, ensure the donor's safety and comfortable.

Ultrafiber leather color table



* Note: the above is the physical scan, there may be a certain color difference in kind.

SKE-136

FEATURES AND HIGHLIGHTS

Technical Parameters

▪ Seat surface size (length×width)	2020mm×600mm(±10mm)
▪ Headrest size	470×260×100mm(±10mm)
▪ Backrest size	610×600mm (±10mm)
▪ Seat size	500×600mm (±10mm)
▪ Leg size	570×600mm (±10mm)
▪ Footrest size	260×640mm (±10mm)
▪ Armrest size	500 ×170mm (±10mm)
▪ Mattress thickness	90mm (±10mm)

Technical Configuration

☑ Hand Controller	1 Pc
☑ Manual Cpr Device	1 Pc
☑ Usb Port	1 Pc
☑ 3 - Inch Silent Brake Wheels	4 Pcs
☑ Side Guide Rails	2 Pcs
☑ Dining Table	1 Pc
☑ External Power Supply	1 Pc
☑ Armrests	2 Pcs
☑ Plug - In Headrest	1 Pc(with a small pillow)
☑ Seamless Soft - Padded Seat Cushion	



Functions

▪ Backrest Lifting	0° - 70° (±5°)
▪ Leg Lifting	- 50° - 0° (±5°)
▪ Height Adjustment	620 - 870Mm (±10Mm)
▪ Armrest Rotation Angle	0° - 270°, Can Be Flipped Backward By 180°
▪ Armrest Height Adjustment	90Mm (±10Mm)
▪ Footrest Adjustment Size	150Mm (±10Mm)
▪ With Manual Cpr Function	



Plastic
Made of brand - new environmentally friendly plastic, non - toxic and odorless. It can be sterilized at temperatures above 100°C with a tensile strength of 30 MPa.



Mechanical function
Includes 5 functions - backrest elevation, leg rest elevation, height adjustment, footrest, and overall reclining up to 6°.




CPR
Quick-release mechanism for CPR, marked with a conspicuous red identifier.



Backrest and leg adjustment
Offers extensive adjustability for backrest and leg rest, allowing easy transition between sitting and lying positions.




Lowest position
Seat can be lowered to 620mm for easy access.



Foot padel adjustment
Removable footrest adjustable forward and backward by 150mm.



Multi-directional armrest
Adjustable, rotatable, and can be flipped up to 180°, allowing patients to find their most comfortable position.



Headrest
Adjustable headrest with 4 locking positions and an adjustment range of 180mm.



Material
Featuring a five-section (headboard, backrest, seat, leg rest, footrest) integrated electric dialysis chair design that conforms to ergonomics, providing comfort for patients. The design includes European prototype elements.



Steels
Made from raw materials of Shanghai Baosteel, a Fortune 500 company.



Seat surface
Uses nano stain-resistant leather with seamless craftsmanship, soft to touch, easy to clean, alcohol-resistant.



Tray table
Equipped with anti-slip edges, cup holder, and phone slot, made of ABS material. Dimensions: 440*330mm.

SKE-188A

ELECTRIC DIALYSIS CHAIR

Technical parameters

▪ Size(LxW)	2040 x 700 mm(±10mm)
▪ Back	820 x 700 mm(±10mm)
▪ Seat	470 x 700 mm(±10mm)
▪ Leg	660 x 700 mm(±10mm)
▪ Foot board	680 x 500 mm(±10mm)
▪ Pillow	360 x 230 mm(±10mm)
▪ Mattress thickness	80 mm(±10mm)











Functions

☒ Back-rest Adjustment	0° - 70°(±5°);
☒ Knee-rest Adjustment	-40° - 0°(±5°);
☒ Height adjustment	580 - 670mm(±10mm);
☒ Foot board adjustable distance	150 mm(±10mm);


Technical configuration

☐ Hand Controller	1pc
☐ Power Line	1pc
☐ Foldable aluminum guardrail	1 set
☐ 4 Silent Castor	4pc
☐ Pillow	1pc





CONSTANT INNOVATION




Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

Portugal Showroom

Av. da Liberdade, nº 101
8150-101
São Brás de Alportel
Portugal
TEL: +351 289 090 212

Brazil Showroom

Rua Lauro Muller 116, Sala
1503 - Botafogo, Rio de
Janeiro - RJ, 22290-160
E-mail: mehecobrasil@
gmail.com

Indian Showroom

1460 M.I.E, Bahadurgarh,
Haryana-124507, INDIA
TEL: +91 96677 70747
E-mail: ind@saikangmedical.
com

China Headquarter

No.35 Lehong Road,
Zhangjiagang City,
Jiangsu Province, China.
TEL: +86 18021231901
E-mail: export@
saikangmedical.com

South East Asia Branch

**SAIKANG Medical
(Thailand) Co., Ltd.**
Ample Tower 909 Room 11/3
Debaratana Road,
Bangna Nuea, Bangkok
10260 Thailand
TEL: +66-936-569469

US Branch

**Healthward Technology,
INC**
17800 castleton st ste 180,
Los Angeles
TEL: +1 (380) 254 4913

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🕒 24/7
✉️ export@saikangmedical.com
🌐 www.saikangmedical.com



Company video



General catalogue