



AIR MATTRESS

SAIKANG MEDICAL

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

Addr. No.35 Lehong Road,Zhangjiagang City,Jiangsu Province,China.

Tel. 0086-18021231901

Fax. 0086-512-58525850

Web. <http://www.jssaikang.com>
<http://www.chinasaikang.com>

E-mail sales@chinasaikang.com




ABOUT US

Jiangsu Saikang Medical Equipment Co.,Ltd. has been established since 2002. It is located in Zhangjiagang City, along the Economic circle of Yangtse River Delta, enjoying the highly developed economic expansion and convenient transportation access.

Our factory is combined with manufacture, sales, service, and supply of Hospital Furniture, Emergency, Surgical products. We are devoted to providing a one-stop sales model to customer; We give more important care on our customer's satisfaction.

Workshop

So far, over 90 countries and regions have enjoyed our products and service. Our success has been mainly achieved by following the corporate philosophy of being firmly committed to "Integrity is the base, Customer is the first".



Bubble Air Mattress



Two alternative colors



Folded



Junction

It can provide excellent pressure reduction through continuous movement. It has been proven clinically beneficial and suitable for all healthcare environment.

Mattress

- Arrangement: 2.5"bubble pad
- Material: medical grade PVC
- Dimension: 190*90*6.4cm
- Weight: 2.6kg

Model: SKP006

Introduction

This therapy system is designed for those with superficial established pressure sore or with existing tissue damaged, or assessed to be a low risk or be at medium risk, even at very high risk of pressure sore development.



Repair Kit




Taiwan NanYa PVC K80



Pump

- Electrical: AC220-240V/100-120V, 50HZ/60HZ
- Air output: $\geq$ 4L/min
- Pressure Range: $\geq$ 12kPa(90mmHg)
- Cycle Time: 5~6min



Model: SKP012



CPR



Cell

Advantages

- One time forming mattress
- Taiwan Nanya PVC,More durable and high elastic
- End flaps 38cm length.
- 100% testing 24 hours before shipment.

Mattress


- Power supply: 110V/220V/50-60Hz
- Air output: 6-8 liter/min.
- Pressure range:50mmHg - 120 mmHg
- Cycle time: 12 min.
- Power consumption: 7W
- Fuse: 1A
- Size: 27.5 (L) x13.5 (W) x 9 (H) cm
- Weight: 1.6kg
- Size: 27.5 (L) x13.5 (W) x 9 (H) cm
- Weight: 1.6kg



Repair Kit

Pump

- Standard Size: 80 "(L) x32-36 "(W)
- 4" 24 cells , 4.5" 18 cells
- Top cover: Nylon PU /strength
- Air cells:Nylon / PVC/TPU
- Bottom: tarpaulin /nylon PU
- Mattress weight : 5.5-6.5kg
- Weight Capacity: 135-145kg
- CPC, Filter, fuse, hook, ventilate option;Aluminum pump system/ AB,ABC alternating option



Taiwan NanYa PVC K80

Tubular Air Mattress

Model: SKP013

Mattress

- Standard Size: 80 "(L) x32-36 "(W)
- 4" 24 cells , 4.5" 18 cells
- Top cover: Nylon PU /strength
- Air cells:Nylon / PVC/TPU
- Bottom: tarpaulin /nylon PU
- Mattress weight : 5.5-7kg
- Weight Capacity: 135-145kg
- CPC, Filter, fuse, hook, ventilate option; Aluminum pump system / AB,ABC alternating option



CPR for quick deflation



Junction and zipper



PU Cover

Pump

- Power supply: 110V/220V/50-60Hz
- Air output: 6-8 liter/min.
- Pressure range:50mmHg - 120 mmHg
- Cycle time: 12 min.
- Power consumption: 7 W
- Fuse: 1A
- Size: 27.5 (L) x13.5 (W) x 9 (H) cm
- Weight: 1.6kg

Foam Base Air Mattress



Model: SKP014

Pump

- Power supply: 110V/220V/50-60Hz
- Air output: 6-8 liter/min.
- Pressure range:50mmHg - 120 mmHg
- Cycle time: 12 min.
- Power consumption: 7 W
- Fuse: 1A
- Size: 27.5 (L) x13.5 (W) x 9 (H) cm
- Weight: 1.6kg



CPR for quick deflation



Junction and zipper

Mattress

- Standard Size: 80 "(L) x32-36 "(W)
 - 5+3" foam pocket
 - Top cover: Nylon PU /strength
 - Air cells:Nylon / PVC/TPU
 - Bottom: tarpaulin /nylon PU
 - Mattress weight : 10kg
 - Weight Capacity: 165kg
- CPC, Filter, fuse, hook, pillow,CPR, ventilate option, zipper; Aluminum pump system / AB,ABC alternating option; High density foam insert foam pocket



PU Cover




Taiwan NanYa PVC


Cell on/in Cell Air Mattress

Pump


- Power supply: 110V/220V/50-60Hz
- Air output: 6-8 liter/min.
- Pressure range:50mmHg - 120 mmHg
- Cycle time: 12 min.
- Power consumption: 7 W
- Fuse: 1A
- Size: 27.5 (L) x13.5 (W) x 9 (H) cm
- Weight: 1.6kg



Taiwan NanYa PVC



Cell on/in Cell



Model: SKP015




Zipper and PU cover



CPR for quick deflation



Junction and zipper



Cells and Pipe

Mattress

- Standard Size: 80 "(L) x32-36 "(W)
 - 8" cell on cell 18/ 20cells; 200*90*20.5cm
 - Top cover: Nylon PU /strength
 - Air cells:Nylon / PVC/TPU
 - Bottom: tarpaulin /nylon PU
 - Mattress weight : 10kg
 - Weight Capacity: 180kg
- CPC, Filter, fuse, hook, pillow, CPR, ventilate option, zipper; Aluminum pump system / AB,ABC alternating option; 8" cell on cell