

24hrs: 0086-18021231901 
www.saikangmedical.com

SAIKANG (Portugal Showroom)
Av. da Liberdade, nº 101 8150-101
São Brás de Alportel Portugal
TEL: +351 289 090 212

SAIKANG (Brazil Showroom)
Rua Lauro Muller 116, Sala 1503 - Botafogo,
Rio de Janeiro - RJ, 22290-160
E-mail: mehecobrasil@gmail.com

SAIKANG (China Headquarter)
No.35 Lehong Road, Zhangjiagang City,
Jiangsu Province, China.
TEL: +86 18021231901
E-mail: export@saikangmedical.com

SAIKANG (South East Asia Branch)
SAIKANG Medical (Thailand) Co., Ltd.
Ample Tower 909 Room 11/3 Debaratana Road,
Bangna Nuea, Bangkok 10260 Thailand
TEL: +66-936-569469

SAIKANG (US Branch)
Healthward Technology, INC
17800 castleton st ste 180, Los Angeles
TEL: +1 (603) 416-3063

JIANGSU SAIKANG MEDICAL EQUIPMENT CO., LTD.
www.saikangmedical.com

Stock Code: 870098

SAIKANG[®]

PARTS FOR HOSPITAL
FURNITURE



Factory Tour

Website

H26-1

LEADING BRAND OF CHINA
MEDICAL FURNITURE




LEADING BRAND OF MEDICAL FURNITURE

Patients in medical rehabilitation often have to face huge challenges. Many people who suffering from diseases or accidents are faced with retraining their minds and bodies and must adapt to the new environment. In this case, the positive atmosphere would offer decisive help. Hospital beds and furniture should provide the best support for patients and nurses. Considering this case, Saikang devoted to provide medical care solutions and perfect choices for medical institutions and rehabilitation centers, to help patients achieve rapid or long-term recovery.

Saikang was founded in 2002, now is a leading global provider of medical furniture and solutions. Head office locates in Suzhou, and branches in US, Dubai, Russia, Brazil, Thailand and Kenya, having 400 employees worldwide. Saikang has formed a large global network of research and development, marketing and service. Currently, Saikang products and medical solutions have been used in more than 149 countries and regions in the world.

V

Series



V21-2
Plastic Head Foot Board
Size: 990×460mm Hook Distance: 740mm




V21-3
Plastic Side Rail
Size: 942/905×350



V21-1
Plastic Bed Platform
Size: 900×1925mm



C5E
5'' Double Side Castors, Central Locking




E8D
Motor
Customized By Specification



Pa4
Four-Hook IV Pole
Size: Φ 16-Φ 19mm / 890-1540mm



Y Series



Y19-1
Plastic Head Foot Board
Size: 920×465mm




Y19-2
Plastic Side Rail
Size: 520/1080×350mm



S4C
PP Bed Platform (Removable and Washable)
Size: 900×1925mm



C5E
5'' Double Side Castors, Central Locking




E8D
Motor
Customized By Specification



Pa4
Four-Hook IV Pole
Size: Φ 16-Φ 19mm / 890-1540mm



D Series



D17-2
Plastic Head Foot Board
Size: 1000×570mm Hook Distance: 760mm




D17-1
Plastic Side Rail
Size: 650/975×340mm



S1F
Powder Coated Steel Bed Platform
Size: 830/900×1925mm



C5E
5'' Double Side Castors, Central Locking




E8D
Motor
Customized By Specification



Pa4
Four-Hook IV Pole
Size: Φ 16-Φ 19mm / 890-1540mm



B Series



102A
Plastic Head Foot Board
Size: 920×560/485mm Hook Distance: 760mm




G3A
Plastic Side Rail
Size: 716×345mm



S1F
Powder Coated Steel Bed Platform
Size: 830/900×1925mm



C5E
5'' Double Side Castors, Central Locking




E8D
Motor
Customized By Specification



Pa4
Four-Hook IV Pole
Size: Φ 16-Φ 19mm / 890-1540mm



K Series



K23-1
Plastic Head Foot Board
Size: 980×480mm Hook Distance: 910mm



G5C
Collapsible Aluminium Alloy Side Rail (6)
Size: 1465×425mm



S1F-1
Powder Coated Steel Bed Platform
Size: 830/900×1925mm



C3A
5'' Covered Castors (Slider)



PF
Folding Crank
Customized By Specification



Pa4
Four-Hook IV Pole
Size: Φ 16-Φ 19mm / 890-1540mm



A

Series



SKB041-1



103A

Plastic Head Foot Board

Size: 950×530/440mm Hook Distance: 910mm



G5C

Collapsible Aluminium Alloy Side Rail (6)

Size: 1465×425mm



S1F

Powder Coated Steel Bed Platform

Size: 830/900×1925mm



SKB041-P5

Bed Plate

Size: 1840×550mm Thickness: 20mm



G4B


Plastic Side Rail

Size: 1070×300mm



C6E

6'' Double Side Castors, Central Locking



C3A

5'' Covered Castors (Slider)



PF

Folding Crank

Customized By Specification



P2L

Dinning Table Top

Size: 815×340mm



M1E

Metal Crank-Rocker

Customized By Specification



Pa4

Four-Hook IV Pole

Size: Φ 16-Φ 19mm / 890-1540mm



SKB041-1-P24

Polyurethane Foam with High Density Sponge

Density: 28kg/m3



D17-2

Plastic Head Foot Board

Size: 1000×570mm Hook Distance: 760mm



102A

Plastic Head Foot Board

Size: 920×560/485mm Hook Distance: 760mm



103A

Plastic Head Foot Board

Size: 950×530/440mm Hook Distance: 910mm



SK-201-C

Plastic Head Foot Board

Size: 920×560mm Hook Distance: 580mm



SK-867

Plastic Head Foot Board

Size: 940×610mm Hook Distance: 760mm



SK-887

Plastic Head Foot Board

Size: 946×520mm Hook Distance: 570mm



Y19-1

Plastic Head Foot Board

Size: 920×465mm



V21-2

Plastic Head Foot Board

Size: 990×460mm Hook Distance: 740mm



K23-1

Plastic Head Foot Board

Size: 980×480mm Hook Distance: 910mm



105A

Plastic and Metal Head Foot Board

Size: 918×500mm Hook Distance: 910mm



T847

Aluminium Alloy Head Foot Board

Size: 910×570/420mm Hook Distance: 830mm



T847-1

Aluminium Alloy Head Foot Board

Size: 760×590/430mm Hook Distance: 640mm



SK-B013

Plastic Head Foot Board

Size: 950×530mm Hook Distance: 900mm



SK-B018

Plastic Head Foot Board

Size: 980×530/440mm Hook Distance: 910mm



102H

Plastic Head Foot Board

Size: 940×545/490mm Hook Distance: 910mm



Q22-2

Plastic Head Foot Board

Size: 843×560mm Hook Distance: 800mm



SK-22

Plastic Head Foot Board

Size: 900×600/500mm Hook Distance: 760mm



SK-876

Plastic Head Foot Board

Size: 950×430/420mm Hook Distance: 760mm



101A

Plastic Head Foot Board

Size: 890×560mm Hook Distance: 470mm



SK-868

Plastic Head Foot Board

Size: 930×560mm Hook Distance: 620mm



SK-885

Plastic Head Foot Board

Size: 980×560mm Hook Distance: 620mm



115A

Wooden Head Foot Board

Size: 980×600/500mm Hook Distance: 830mm



115B

Wooden Head Foot Board

Size: 1100×600mm Hook Distance: 1050mm



104D

Plastic and MDF Head Foot Board

Size: 1000×500/430mm Hook Distance: 820mm



D17-1
Plastic Side Rail
Size: 650/975×340mm



G2B
Plastic Side Rail
Size: 840/1100×310mm



G3A
Plastic Side Rail
Size: 716×345mm



Y19-2
Plastic Side Rail
Size: 520/1080×350mm



T20
Plastic Side Rail
Size: 1093/820×350mm



V21-3
Plastic Side Rail
Size: 942/905×350



SK-863
Plastic Side Rail
Size: 800×365mm



SK-865
Plastic Side Rail
Size: 560/1050×380mm



SK-892
Plastic Side Rail
Size: 1030/860×380



SK-862
Plastic Side Rail
Size: 740×340mm



G1C
Plastic Side Rail
Size: 720×300mm



G4B
Plastic Side Rail
Size: 1070×300mm



SK-861
Plastic Side Rail
Size: 710×350mm



SK-402
Aluminium Alloy Tuck Away Side Rail
Size: 775×570mm



G8C
Wooden Side Rail
Size: 1920×100mm



G6D
Collapsible Bending Tube Side Rail (4)
Size: 950×430mm



G6D
Collapsible Bending Tube Side Rail (5)
Size: 1205×430mm



G6D
Collapsible Bending Tube Side Rail (6)
Size: 1465×430mm



G5C
Collapsible Aluminium Alloy Side Rail (4)
Size: 950×425mm



G5C
Collapsible Aluminium Alloy Side Rail (5)
Size: 1205×425mm



G5C
Collapsible Aluminium Alloy Side Rail (6)
Size: 1465×425mm



G5C-1
Collapsible Aluminium Alloy Side Rail (4)
Size: 960×425mm



G5C-1
Collapsible Aluminium Alloy Side Rail (5)
Size: 1220×425mm



G5C-1
Collapsible Aluminium Alloy Side Rail (6)
Size: 1480×425mm



G6C
Collapsible Stainless Steel Side Rail (4)
Size: 930×450mm



G6C
Collapsible Stainless Steel Side Rail (5)
Size: 1200×450mm



G6C
Collapsible Stainless Steel Side Rail (6)
Size: 1470×450mm




G5C
Collapsible Painted Wooden Color Side Rail
Size: 1465×425mm




G12C
Collapsible Full-length Side Rail
Size: 1830×425mm




SK-B12
Three Sections Collapsible Side Rail
Size: 1360×490mm




C6E
6" Double Side Castors, Central Locking




C5B
5" Hard Shell Castors




C3A
5" Covered Castors (Slider)




C3A
4" Covered Castors




C4B
3" Noiseless Castors




C5E
5" Double Side Castors, Central Locking




C4B
5" Noiseless Castors




C3A
5" Covered Castors (Screw)




C3A
3" Covered Castors




C5F
2" Noiseless Castors




C2F
5" Single Side Castors, Central Locking




C5D
5" Double Side Castors




C3A
4" Covered Castors (Screw)



C3F
3" Noiseless Castors



C5G
2" Noiseless Castors




S1F
Powder Coated Steel Bed Platform
Size: 830/900×1925mm



S4A
Phenolic Resin Bed Platform (X-ray Translucent)
Size: 830/900×1925mm



S4D
Plastic Bed Platform
Size: 900×1895mm




P2L
Dinning Table Top
Size: 815×340mm




SKH046-SC
Telescopic Dinning Table
Size: (890-1130)×350mm



S1G
10 Pieces Powder Coated Steel Bed Platform
Size: 827/895×1993mm



S4B
Powder Coated Steel Bed Platform
Size: 830/900×1925mm



V21-1
Plastic Bed Platform
Size: 900×1925mm



P2R
Moving Dinning Table Top
Size: 890×400mm



S1A
Plastic Bed Platform
Size: 830/900×1935mm



S4C
PP Bed Platform (Removable and Washable)
Size: 900×1925mm



Q22-1
Plastic Bed Platform (For Children bed)
Size: 780×1730mm



P2S
Moving Dinning Table Top
Size: 890×440mm



M1E
Metal Crank-Rocker
Customized By Specification



PF
Folding Crank
Customized By Specification



E8D
Motor
Customized By Specification



E8D
Motor
Customized By Specification



M1D
Metal Crank Handle
Size: 220×150mm



M1F
Crank Handle
Size: 165×250mm



E8D
Motor
Customized By Specification



E8F
Control Box
Customized By Specification



E8F
Control Box
Customized By Specification



E8G
Back up Battery
Customized By Specification



E5a
Portable Control Board
Customized By Specification



E5a
Hand Controller
Customized By Specification



E5a
Hand Controller
Customized By Specification



Pa2
Four-Hook IV Pole
Size: Φ 16-Φ 19mm / 910-1500mm



P3C
Orthopaedics Balkan Frame
Size: 2100×1500×1000mm



P3F
Lifting Handle
Size: 100×100×650mm



P5F
Medical Record Card Socket
Size: 137×75mm



SKR023
Medical Record Clip
Size: 320×230mm



Pa2
Two-Hook IV Pole
Size: Φ 16-Φ 19mm / 910-1500mm



P3C-1
Double-arm Lifting Pole
Size: 900×1500mm



P3E
Lifting Handle
Size: 250×145×590mm



P5G
Writing Board
Size: 300×250×18mm



P5H
Bumper Wheel
Size: Φ 100mm



Pa1
IV Pole
Size: Φ 16-Φ 19mm / 910-1500mm



P3A
Single-arm Lifting Pole
Size: 1500mm




P5D
48 Lattice Medicine Box
Size: 510×360×60mm



X01-101
Baby Basin
Size: 805×455×255mm



P5K
Bumper Strip
Size: 30×5mm



DESIGN & MOLDING SERVICE PROCESS

01

Customer Concept Sketch


- * Free of charge
- * Initial idea or hand-drawn sketch from customer



03

Detailed Design Drawing

- * Requires 30% advance payment
- * Refined and accurate design draft (10 days)



05

Mold Fabrication


- * Includes mold trial and up to three free modifications
- * Requires 40% final payment
- * Lead time: 45 days



02

Preliminary Rendering

- * Free of charge
- * Rough simulation drawing (5 days)



04

3D Drawing Confirmation (PPT Presentation)

- * Requires additional 30% payment
- * 3D design presented in PPT format for customer review and confirmation
- * After confirmation, 3D prototype sample will be produced for reporting and validation (15 days)



06

First Batch Sample Production

- * Mass production of the first sample batch after mold confirmation